

FREE PROGRAMMABLE INPUTS:


PROGR. INPUT

(PIN3 (+) / PIN 25) but is not wired yet!

ADDITIONAL INPUT 1

ADDITIONAL INPUT 2

ADDITIONAL INPUT 3



There is **1** programmable input in the basic version.

With the option "I/O extension board" further **3** inputs are available.

Each input can be programmed individually according to the following list:

- **NO FUNCTION**

reset of setting

- **PROGRAM START PULSE**

if selected, program x can be started (**external start**)

- **PROGRAM STOP PULSE**

if selected, program x will be stopped (terminate with ESC)

- **PROGRAM RUN DURING PULSE**

A program is executed as long as there is a continuous signal. If the signal drops, the program will be stopped.

- **BOTTLE CHANGE PULSE**

Pulse signal: **<= 3sec** means "advance to next bottle"

>= 5 sec means "move to bottle No. 1"

- **SAMPLING PULSE**

A pulse triggers the sample extraction

The pulse signal has to be **> 50 ms!**

Remark: This function is only possible if **no** program is running. In this case the unit is controlled externally (e.g. via SPS).

- **MANUAL SAMPLE**

A manual sample is triggered.

There is no registration in the info memory (ideal for official samples or tests).

- **ACKNOWLEDGE ERROR**

Accumulated error messages can be acknowledged.